



Orientation Discrimination: Lesson 3

In this lesson, you will learn how to:

- create noise masks,
- use a timeline event to specify the temporal relationship between the stimuli events.

This lesson assumes you have been through [Lesson 1](#) of this tutorial and became familiar with adding and inspecting the various types of events.

Difficulty: 1/5

Duration: 15 mn

TEMPORAL MASKING

This lesson makes the basic orientation discrimination task more elaborate by displaying 2D noise masks before and after the target stimulus (aka sandwich paradigm with forward and feedback masking stimuli).

First, duplicate the original experiment and rename the copy **Orientation Discrimination 3** as you did in [Orientation Discrimination: Lesson 2](#).

Step 1: Adding a Timeline Event

Option-click on the arrow in front of the **Orientation Discrimination 3** experiment to reveal its whole hierarchy.

Title	Category		
Orientation Discrimination 3	Experiment	⌵	🔒
Staircase	Method	⌵	🔒
2IFC	Procedure	⌵	🔒
Gabor	Visual Stimulus	⌵	🔒
Orientation Discrimination 2	Experiment	⌵	🔒
Orientation Discrimination	Experiment	⌵	🔒
Storage		⌵	
Trash		⌵	

Select the **2IFC** procedure and insert a new group event using the '+' folder icon. Change the sub-category and title of this new group event to **Composed Stimuli/Timeline** and **Timeline**, respectively.

Title	Category		
Orientation Discrimination 3	Experiment	⌵	🔒
Staircase	Method	⌵	🔒
2IFC	Procedure	⌵	🔒
Timeline	Composed Stimuli	⌵	🔒
Gabor	Visual Stimulus	⌵	🔒
Orientation Discrimination 2	Experiment	⌵	🔒
Orientation Discrimination	Experiment	⌵	🔒
Storage		⌵	
Trash		⌵	



Move the **Gabor** stimuli **INSIDE** the **Timeline** event so it gets indented to the right as illustrated.

Title	Category	Settings	Lock
Orientation Discrimination 3	Experiment	⌵	🔒
Staircase	Method	⌵	🔒
2IFC	Procedure	⌵	🔒
Timeline	Composed Stimuli	⌵	🔒
Gabor	Visual Stimulus	⌵	🔒
Orientation Discrimination 2	Experiment	⌵	🔒
Orientation Discrimination	Experiment	⌵	🔒
Storage		⌵	
Trash		⌵	

Select the **Timeline** event and click on the **Inspect** button to inspect its properties. If you did not select yet a sub-category for this **Composed Stimuli** event above, click on the **Timeline** tab now (1).

The **Timeline** properties present a list of embedded stimulus events: so far there is only one, the **Gabor** stimulus.

Click on the **OK** button to close the properties panel (2). The next step is to add masking stimuli before and after the **Gabor** stimulus.

Step 2: Adding Masking Stimuli

Select the **Timeline** event, and insert a new stimulus event using the '+' file icon. Change the sub-category and title of this new group event to **Visual Stimulus/Grating** and **Forward Mask**, respectively.

Title	Category	Settings	Lock
Orientation Discrimination 3	Experiment	⌵	🔒
Staircase	Method	⌵	🔒
2IFC	Procedure	⌵	🔒
Timeline	Composed Stimuli	⌵	🔒
Forward Mask	Visual Stimulus	⌵	🔒
Gabor	Visual Stimulus	⌵	🔒
Orientation Discrimination 2	Experiment	⌵	🔒
Orientation Discrimination	Experiment	⌵	🔒
Storage		⌵	
Trash		⌵	

Tip: You could also duplicate the **Gabor** event, rename the copy **Forward Mask**, and move it directly above the original **Gabor**.



Select the **Forward Mask** event and edit its properties by clicking on the **Inspect** button (or press **Apple-i**).

Select the **Grating** tab if necessary.


Using the pop-up menu, set the carrier type to **Noise**; then select the **2D**, **Gaussian**, and **White** noise options.

Set the envelope to **Hard Edge** with its radius at **2** deg (this must be larger than the sigma used for the **Gabor** stimulus).

Click on the **OK** button to validate the changes.

Duplicate the **Forward Mask** stimulus, rename the copy to **Backward Mask**, and move it to just after the **Gabor** stimulus inside the **Timeline** event, as illustrated.

The next step is to edit the **Timeline** properties to specify the temporal relationship between the masks and the target stimulus.




Title	Category		
ψ Orientation Discrimination 3	Experiment	✓	🔒
Staircase	Method	✓	🔒
2IFC	Procedure	✓	🔒
Timeline	Composed Stimuli	✓	🔒
Forward Mask	Visual Stimulus	✓	🔒
Gabor	Visual Stimulus	✓	🔒
Backward Mask	Visual Stimulus	✓	🔒
ψ Orientation Discrimination 2	Experiment	✓	🔒
ψ Orientation Discrimination	Experiment	✓	🔒
Storage			
Trash			

Step 3: Customizing the Timeline

Select the **Timeline** event, and click on the **Inspect** button to edit its settings.

There are now three stimuli in the **Timeline** table. The green area in the table indicates when and how long each stimulus is displayed. By default and as illustrated, they are presented simultaneously for the same length of time.





Note that the total duration of the timeline, as well as the onset and duration for the selected stimulus, are indicated above the table.

Select the **Gabor** entry in the table and set its **Onset** column to **ISI** using the pop-up menu (**ISI** = inter-stimulus interval; **SOA** = stimulus-onset asynchrony).

#	Events	Onset	To	Timeline
0	Forward	SOA	NA	500 ms
1	Gabor	ISI	NA	500 ms
2	Backward	SOA	NA	500 ms

Click on **NA** (not available) in the **to** column and select the entry **0: Forward Mask** from the pop-up menu. This selection ensures that the **ISI** for the **Gabor** stimulus occurs relative to the **Forward Mask** stimulus.

- ✓ NA
- 0: Forward Mask
- 1: Gabor
- 2: Backward Mask

The **Gabor** timeline is now split in two: a black zone corresponding to the period before the **ISI** is applied, and a green zone during which the stimulus is displayed.

#	Events	Onset	To	Timeline
0	Forward	SOA	NA	500 ms
1	Gabor	ISI	0:...	500 ms
2	Backward	SOA	NA	500 ms

Note that the **ISI** onset is set to 0 second by default which means:

- the **Gabor** onset occurs immediately at the offset of the **Forward Mask**,
- the total duration has also been automatically updated (from 0.5 to 1 second) to account for the extra time introduced by the **Gabor ISI** relative to the **Forward Mask**.



Select the **Backward Mask** entry in the table and set its **Onset** column to **ISI** and its **to** column to **1: Gabor**. This specifies an **ISI** (of 0 second by default) between the **Gabor** stimulus and the **Backward Mask** stimulus.

The **Gabor** stimulus is now embedded in a sandwich masking paradigm.

#	Events	Onset	To	Timeline
0	Forward	SOA	NA	500 ms
1	Gabor	ISI	0:...	500 ms
2	Backward	ISI	1:...	500 ms

Click on the **OK** button to validate the changes and return to the **Designer** panel.
Check and run the experiment!

Conclusion

In this lesson you learned how to create noise masks and use a timeline to precisely specify the relative timing of the stimuli, particularly their presentation order through the use of an ISI.

In the next lessons, you will learn how to:

- add a spatial and temporal context ([Lesson 4](#)),
- add motion & dynamics through the use of a 1st-order drifting Gabor ([Lesson 5](#)).

If you are interested in further exploring masking effects on orientation discrimination, you may try to extend this lesson by investigating the effects of various masking parameters, such as noise contrast, noise size, noise duration or ISI duration. In that case, make sure to go through [Lesson 2](#) of this tutorial as it shows how to configure different experimental conditions and run interleaved staircase methods.